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Destinazione d'uso

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enINSTRUCTIONS FOR USE

Intended use

Masks for CPAP (continuous positive airway pressure) ventilation in a hospital environment.

Indications for use

Indicated for the treatment of acute respiratory insufficiency of: hypoxemic patients (acute cardiogenic pulmonary oedema - in CPAP) and post-operative atelectasis. For other pathologies the treatment may only be performed under careful monitoring by expert operators.

Restrictions of use

- The device must be used by qualified and trained medical/nursing staff.
- To be used for administration of air and oxygen.
- To be used on adult patients weighing > 30 kg.
- If used properly, the device can be used continuously for a maximum of 7 days after which it must be replaced.

Contraindications

- Coma
- uncooperative patient
- cardiac arrest
- hemodynamic instability
- recent oesophageal and gastro-surgical operations
- heavy bleeding of the upper digestive tract
- obstruction of the upper airways
- pneumothorax
- impossibility to remove secretions
- Cranial/facial traumas or burns
- Vomiting
- Not collaborating patients and unable to remove the mask themselves.

Complications/side effects

Facial irritations and traumas. Discomfort during prolonged use. Flow leaks from interface. Eye irritation / conjunctivitis.

Instructions for use

Remove the device from the package. Attach the connector to the ventilation circuit. Check that the connection is secure before proceeding with administration of the therapy gases. The mask may be equipped with a PEEP valve. With adjustable PEEP: set the desired pressure by acting on the cap (clockwise to increase the PEEP, anticlockwise to decrease it). Possible adjustments from 0 - 20 cm H2O (mbar). With pre-calibrated PEEP: Always check that the pressure indicated on the device is as desired. The mask is fitted with an atraumatic cushion: always check inflation before proceeding with gas delivery.

Warnings and precautions

⚠ Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients.

- In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.

- Monitor the patient's clinical parameters during the treatment.
- If the patient's conditions deteriorate or there is no improvement within the scheduled times, it is recommended to evaluate an alternative ventilation technique.
- The efficacy of the therapy is considerably affected by the PEEP values set. Therefore, the medical staff are advised to carefully evaluate the pressure level most suited to the clinical condition of the patient. Using too low pressure may not be sufficient for alveolar recruitment. Using too high pressure may be a source of alveolar over distension.
- When used with humidified ventilation gas, check that any condensate that forms does not discomfort the patient.
- The patient must REMOVE all accessories such as glasses, earrings, hair clips, combs and any other metal objects.
- The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size.
- Use only on unbroken skin.
- The device can be used for maximum 7 days. Depending on the clinical picture of the patient the medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.
- Non-sterile. Do not sterilize.
- Do not use on another patient.
- Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°C/+50°C).
- Dispose of the materials immediately after use conforming to the current laws and regulations.
- Use only with devices that have connections in compliance with the ISO 5356-1 standards.
- Non invasive ventilation shall be interrupted if there is one of the following conditions:
 - worsening state of consciousness and respiratory distress
 - loss of airway protection
 - unchanged PaCO₂ (during two subsequent EGA performed after max 1 hour)
 - persistence of severe hypoxemia
 - severe and uncontrollable hemodynamic instability
 - patient / ventilator synchronization problems
 - uncontrollable secretions
 - interface intolerance

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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