en INSTRUCTIONS FOR USE

The instruction for use (IFU) is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE Distribute instructions to all product locations and users.

These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using the product. Failure to read the instructions, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury. **CAUTION** A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

INTENDED USE To make secure, gas tight connections to other respiratory standard taper connectors and/or provide a respiratory pathway between a breathing system and patient's airway, facemask, face mask, monitoring ports.

INTENDED USER For use by or under the supervision of appropriately qualified clinician or healthcare professional in a hospital and homecare environment.

LIMITATIONS OF USE WARNING Do not use with equipment that does not comply to international standards.

PATIENT CATEGORY The product range is suitable for use with all patient populations including adults, paediatric and neonates.

PRE-USE CHECK

1. Prior to attachment, check that all breathing system components are free of any obstructions or foreign bodies and that there is no patient air leak.

2. Connect/disconnect the tapered connections of the breathing system components using a push and twist action (fig. 1).

3. Following attachment the breathing system and all accessories must be checked for function, leakage and occlusions prior to use and that all connections are secure.

IN-USE CHECK

1. Monitor the clinical parameters of the patient during the treatment.

2. Monitor the breathing system and related components throughout use and replace if visual contamination or increased resistance occurs.

WARNING

1. Prior to use consult the respective documentation and instructions for use of all connected devices, or combination of devices and perform system tests as required.

2. Single patient use product. Do not reuse. Reuse poses a threat to patient safety in terms of cross infection and device malfunction and breakage. Cleaning can lead to device malfunction.

3. Do not use product past the expiry date.

CAUTION

1. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged or tampered with.

2. The products have been designed, validated and manufactured for single use and for a period not exceeding 7 days.

3. The product must not be re-used, cleaned or sterilized.

4. FOR USA: Rx ONLY: Federal law restricts this device to sale by or on the order of a physician.

STORAGE CONDITIONS Recommended storage at room temperature over the indicated product shelf-life.

DISPOSAL Following use, the product must be disposed of in accordance with local health, infection control and waste disposal regulations.

APPROPRIATE STANDARDS This product complies with relevant requirements of the following standards:

- BS EN ISO 5367: Anaesthetic and respiratory equipment. Breathing tube and sets
- BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets

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