
Bilevel breathing systems – Smoothbore 22mm systems


en	Bilevel breathing systems – Smoothbore 22mm systems
fr	Circuits respiratoires VDNP – Circuits Intérieur lisse 22mm
de	Bilevel Beatmungssysteme – Smoothbore 22mm Systeme
es	Circuitos respiratorios Binivel – circuitos Smoothbore de alma lisa de 22mm
pt	Sistemas respiratórios de dois níveis – Sistemas Smoothbore de 22 mm
it	Circuiti per ventilazione bilevel – Circuiti Smoothbore da 22mm
nl	Bilevel beademingssystemen - 22mm Smoothbore systemen
no	Bilevel pustesystemer – Smoothbore 22mm systemer
fi	Bilevel-hengitysjärjestelmät - Smoothbore 22mm -järjestelmät
sv	Bilevel-andningssystem – Smoothbore 22mm-system
da	Bilevel respirationssystemer – Smoothbore sæt 22mm
el	Bilevel αναπνευστικά συστήματα – Smoothbore συστήματα 22mm
lt	Bilevel kvėpavimo sistemas- 22mm Smoothbore sistemas
pl	Systemy oddechowę typu Bilevel - Systemy gładkoramkowe 22 mm
ru	Дыхательные системы для вентиляции с двухуровневым давлением Bilevel- Гладкоствольные дыхательные системы 22ММ
cs	Bilevel dýchací systémy – 22mm systémy Smoothbore s hladkým vnitřním průřezem
hu	Kétszintű lélegeztető rendszerek – Smoothbore 22mm rendszerek
sl	Dvonivojski dihalni sistemi – Smoothbore 22 mm sistemi
lv	Divlūmeņš elpošanas sistēmas- Smoothbore 22mm sistēmas
et	Bilevel hingamissüsteem-Smoothbore 22 mm süsteemid
bg	Дихателни системи на две нива – Smoothbore 22 мм системи
sr	Dvostepeni sistemi za disanje – sistemi glatke cevi od 22 mm
ro	Circuite de ventilație Bilevel – circuite 22mm Smoothbore
sk	Bilevel dýchacie systémy – Smoothbore 22mm systémy
hr	Dvorazinski sustavi za disanje – Smoothbore sustavi od 22 mm
tr	Bilevel solunum sistemleri – Smoothbore 22mm sistemleri

Figure 1. System set-up / Figure 1. Configuration du système / Abbildung 1. Systemaufbau / Figura 1. Configuración del sistema / Figura 1. Configuração do sistema / Figura 1. Configurazione del sistema / Figur 1. Systemet installen / Figur 1. Systemoppsett / Kuva 1. Järjestelmän määrittäminen / Figur 1. Systeminställning / Figur 1. Opsætning af systemet / Εικόνα 1. Ρύθμιση συστήματος / 1 pav. Sistemos nustatymas / Rysunek 1. Konfiguracja systemu / Рисунок 1. Настройка системы / Obr. 1. Sestavení systému / 1. ábra. A rendszer felépítése / 1. attéls. Sistēmas iestatīšana / Joonis 1. Süsteemi seadistamine / Фигура 1. Конфигурация на системата / Slika 1. Postavljanje sistema / Figura 1. Configurarea sistemului / Obrázok 1. Zostavenie systému / Slika 1. Postavljanje sustava / Şekil 1. Sistemin kurulumu

* Bilevel Vent / Ouverture Bi-Level / Bilevel-Beatm. / Ventilación Binivel / Vent. Bilevel / Sfiato Bilevel / Bilevel-Ventilatie / Bilevel Apning / Bilevel-Hengityskone / Ventilator Med Dubbelnivå / To-Vejs Vent / Αναπνευστικός Διο Επιπέδων / Bilevel Vent. / Respirator Do Terapii Dwupoziomowej / Вентилятор Bilevel / Dvojúrovňový Ventilátor / Kétszintű Lélegeztető / Ventilator Bilevel / Bilevel Ventilators / Bilevel Vent / Вентилятор На Две Нива / Respirator Na Dwa Nivoa („Bilevel“) / Ventilator Bilevel / Dvojúrovňový Ventilátor / Respirator Bilevel / İki Seviyeli Solunum Cihazı / Bilevel Vent

Figure 2. Leak port / Figure 2. Port de détection des fuites / Abbildung 2. Leckageport / Figura 2. Puerto de fuga / Figura 2. Válvula de exalação / Figura 2. Condotto diffusore / Figur 2. Lekaansluiting / Figur 2. Lekkasjeport / Figure 2. Lekkasjeport / Kuva 2. Vuotoliitäntä / Figur 2. Lackageport / Figur 2. Leckageport / Εικόνα 2. Ούρα διφρροής / 2 pav. Nuotėkio anga / Rysunek 2. Port przeciekowy / Рисунок 2. Отверстие для стекания / 2. ábra. Szivárgási nyílás / Slika 2. Odrpina za pušcanje / 2. attéls. Norpūdes pieslēgvieta / Joonis 2. Lekkeport / Фигура 2. Порт за изтичане / Slika 2. Port za curenje / Figura 2. Port de scurgere / Obrázok 2. Odtokový port / Slika 2. Otvor za istjecanje / Şekil 2. Sızdırma portu

* Bilevel Vent / Ouverture Bi-Level / Bilevel-Beatm. / Ventilación Binivel / Vent. Bilevel / Sfiato Bilevel / Bilevel-Ventilatie / Bilevel Apning / Bilevel-Hengityskone / Ventilator Med Dubbelnivå / To-Vejs Vent / Αναπνευστικός Διο Επιπέδων / Bilevel Vent. / Respirator Do Terapii Dwupoziomowej / Вентилятор Bilevel / Dvojúrovňový Ventilátor / Kétszintű Lélegeztető / Ventilator Bilevel / Bilevel Ventilators / Bilevel Vent / Вентилятор На Две Нива / Respirator Na Dwa Nivoa („Bilevel“) / Ventilator Bilevel / Dvojúrovňový Ventilátor / Respirator Bilevel / İki Seviyeli Solunum Cihazı / Bilevel Vent

en The instruction for use is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request. A copy of the instructions for use is also available on the website at: <https://customers.intersurgical.com/ProductDocumentation/>.
NOTE Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

DESCRIPTION

Bilevel breathing system.

INTENDED USE

To deliver and remove respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

INTENDED USER

For use by qualified medical personnel. This product may be used in a homecare setting. Appropriately qualified medical personnel to instruct homecare users on correct usage.

INDICATIONS

The breathing system should be connected between the ventilator and the patient airway connection on a single patient up to a maximum of 7 days. The start of use is defined from the first connection to the equipment. For use in a hospital environment. For use in a homecare environment.

CONTRAINDICATIONS

No known contraindications.

ACCESSORIES

Water traps and breathing filters, HMEs, HMEFs (if included in the configuration)

Follow the instructions in individual IFUs supplied.

Monitoring lines (if included in the configuration)

Monitoring lines provide a leak free connection between the breathing system and a respiratory gas monitor via a luer lock fitting. The connection should be checked prior to use. The line should be checked to confirm it is kink free and a patent pathway exists.

PATIENT CATEGORY

Intended for use with adult and paediatric patients.

PREPARATION (see fig. 1)

1. Attach breathing system limb (A) to the ventilator gas outlet.
2. If a filter is required, place the filter (B) between the ventilator gas outlet and the breathing system. If a machine end filter is required please use Intersurgical Flo-Guard, or the filter supplied with this breathing system.
3. Attach pressure monitoring line (C) to the ventilator pressure monitoring port (if applicable).
4. Remove red dust cap (D) from the breathing system and attach the breathing system to patient interface.

PRE-USE CHECKS

WARNING

1. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.
2. Do not use if package is damaged.
3. Every component of the breathing system must be visually inspected and checked for function, leakage and occlusions, immediately before use on each patient. Ensure that connections are secure.
4. Refer to equipment manufacturers guidelines for machine settings and pre-use check requirements.
5. After full assembly of the breathing system and associated accessories perform a self-test of the ventilator, following the equipment manufacturer's instructions, before use on each patient, if required.
6. Ensure the leak port is not occluded (if applicable). See fig 2.
7. There is a risk of fire associated with oxygen equipment and therapy. Do not use near sparks or open flames.
8. Do not use product past the expiry date.

CAUTION

1. When connecting or disconnecting the system, use a push and twist action.
2. If a machine end filter is required please use Intersurgical Flo-Guard, or the filter supplied with this breathing system.
3. Use the red dust cap to occlude the breathing system patient connection port as and when required in the ventilator self-test protocol.

Ensure red dust cap is held securely to prevent unwanted disconnection during self-test.

IN-USE CHECKS

WARNING

1. Do not stretch, occlude, cut or reconfigure the tubing.
2. Monitor the clinical parameters of the patient during the treatment.
3. Ventilator accuracy can be affected by use of a gas driven nebuliser.
4. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the ventilator performance.
5. This system is not suitable for use with flammable anaesthetic agents.
6. Monitor the breathing system throughout use and replace if increased resistance or contamination occurs.
7. Breathing systems, their parts and accessories are validated for use with specific ventilators. Incompatible parts can result in degraded performance. The responsible organisation is responsible for ensuring the compatibility of the ventilator and all of the parts used to connect to the patient before use.
8. Use of this breathing system with equipment not compliant with appropriate International Standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure compatibility.
9. Leakage of high levels of oxygen at the patient interface may cause patient harm.
10. The materials used in the airways and related equipment may not be compatible with anaesthetic or respirable gases, solutions/suspensions/emulsions that have not been evaluated.
11. Monitor breathing system condensation and clear any mobile condensate that may build up in the tubes during treatment. Ensure condensate is drained towards the machine and away from the patient.
12. If a luer port connector is included in the system, it is for the exclusive use of monitoring respiratory gases. If applicable, check that the cap is securely attached to prevent inadvertent leakage when the port is not in use. Do not introduce gases or liquids via the misconnection of the luer lock port as this can result in serious injury to the patient.
13. Consult individual product labelling to determine if device is MR safe or unsafe.
14. Consult individual product labelling for presence or absence of phthalates.

CAUTION

When the breathing system is not in use, the patient connection should be protected using the red dust cap supplied, to prevent contamination.

DURATION OF USE

Single use product. Maximum duration of use 7 days. Do not reuse. Reuse poses a threat to patient safety in terms of cross infection and device malfunction and breakage. Do not clean. Cleaning can lead to device malfunction and breakage.

CAUTION

1. The breathing system is only to be used under medical supervision by suitably or appropriately qualified medical personnel.
2. This product has been verified for a maximum duration of 7 days use.
3. Following use, the product must be disposed of in line with the local healthcare, hygiene and waste disposal protocols and regulations.

NOTE

The length referenced in the product description on the product label denotes the minimum working length of the breathing system. The full breathing system length is shown in the technical data section of the product label.

WARNING

Ensure that the breathing system technical and performance data are appropriate to the intended patient and treatment.

Patient ventilation and vital signs should be monitored. Please refer to product label for full information.

APPLICABLE STANDARDS

This product complies with relevant requirements of the following standards:

- BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets
- BS EN ISO 5367: Anaesthetic and respiratory equipment. Breathing sets and connectors
- BS EN ISO 80601-2-12: MEE Critical Care Ventilators

For further information, contact Intersurgical - info@intersurgical.com

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